
FabIO Documentation

Release 0.1.4

**H. Sorensen, E. Knudsen, J. Wright,
G. Goret, B. Pauw and J. Kieffer**

March 31, 2014

CONTENTS

1	Getting Started	3
1.1	Introduction	3
1.2	FabIO Python module	3
1.3	Usage	5
1.4	Future and perspectives	5
1.5	Conclusion	5
2	Installation	9
2.1	Dependencies	9
2.2	Installation from sources	9
2.3	Development versions	10
2.4	Test suite	10
3	Changelog	13
3.1	From FabIO-0.1.3 to FabIO-0.1.4:	13
3.2	From FabIO-0.1.2 to FabIO-0.1.3:	13
3.3	From FabIO-0.1.1 to FabIO-0.1.2:	13
3.4	From FabIO-0.1.0 to FabIO-0.1.1:	14
3.5	From FabIO-0.0.8 to FabIO-0.1.0:	14
3.6	From FabIO-0.0.7 to FabIO-0.0.8:	14
3.7	From FabIO-0.0.6 to FabIO-0.0.7:	14
3.8	From FabIO-0.0.4 to FabIO-0.0.6:	14
4	FabIO Package	15
4.1	fabio Package	15
4.2	fabio.fabioimage Module	15
4.3	fabio.fabioutils Module	17
4.4	fabio.file_series Module	19
4.5	fabio.openimage Module	22
4.6	fabio.adscimage Module	22
4.7	fabio.binaryimage Module	22
4.8	fabio.bruker100image Module	23
4.9	fabio.brukerimage Module	23
4.10	fabio.cbfimage Module	24
4.11	fabio.dm3image Module	25
4.12	fabio.edfimage Module	26
4.13	fabio.fit2dmaskimage Module	29
4.14	fabio.fit2dspreadsheetimage Module	29
4.15	fabio.GEimage Module	29

4.16	<code>fabio.hdf5image</code> Module	30
4.17	<code>fabio.HiPiCimage</code> Module	31
4.18	<code>fabio.kcdimage</code> Module	31
4.19	<code>fabio.mar345image</code> Module	31
4.20	<code>fabio.mrcimage</code> Module	32
4.21	<code>fabio.marccdimage</code> Module	32
4.22	<code>fabio.OXDimage</code> Module	32
4.23	<code>fabio.pilatusimage</code> Module	33
4.24	<code>fabio.pixiimage</code> Module	33
4.25	<code>fabio.pnmimage</code> Module	34
4.26	<code>fabio.raxisimage</code> Module	34
4.27	<code>fabio.tifimage</code> Module	35
4.28	<code>fabio.xsdimage</code> Module	35
4.29	<code>fabio.compression</code> Module	36
4.30	<code>fabio.converters</code> Module	38
4.31	<code>fabio.datIO</code> Module	38
4.32	<code>fabio.TiffIO</code> Module	38
4.33	<code>fabio.readbytestream</code> Module	39
5	Indices and tables	41
	Python Module Index	43
	Index	45

Contents:

GETTING STARTED

FabIO is a Python module for reading and handling data from two-dimensional X-ray detectors.

FabIO is a Python module written for easy and transparent reading of raw two-dimensional data from various X-ray detectors. The module provides a function for reading any image and returning a `fabioimage` object which contains both metadata (header information) and the raw data. All `fabioimage` object offer additional methods to extract information about the image and to open other detector images from the same data series.

1.1 Introduction

One obstacle when writing software to analyse data collected from a two-dimensional detector is to read the raw data into the program, not least because the data can be stored in many different formats depending on the instrument used. To overcome this problem we decided to develop a general module, FabIO (FABle I/O), to handle reading and writing of two-dimensional data. The code-base was initiated by merging parts of our `fabian` imageviewer and `ImageD11` peak-search programs and has been developed since 2007 as part of the TotalCryst program suite for analysis of 3DXRD microscopy data. During integration into a range of scientific programs like the FABLE graphical interface, EDNA and the fast azimuthal integration library, `pyFAI`; FabIO has gained several features like handling multi-frame image formats as well as writing many of the file formats.

1.2 FabIO Python module

Python is a scripting language that is very popular among scientists and which also allows well structured applications and libraries to be developed.

1.2.1 Philosophy

The intention behind this development was to create a Python module which would enable easy reading of 2D data images, from any detector without having to worry about the file format. Therefore FabIO just needs a file name to open a file and it determines the file format automatically and deals with `gzip` and `bzip2` compression transparently. Opening a file returns an object which stores the image in memory as a 2D NumPy array and the metadata, called header, in a Python dictionary. Beside the data and header attributes, some methods are provided for reading the previous or next image in a series of images as well as jumping to a specific file number. For the user, these auxiliary methods are intended to be independent of the image format (as far as is reasonably possible).

FabIO is written in an object-oriented style (with classes) but aims at being used in a scripting environment: special care has been taken to ensure the library remains easy to use. Therefore no knowledge of object-oriented programming is required to get full benefits of the library. As the development is done in a collaborative and decentralized way; a comprehensive test suite has been added to reduce the number of regressions when new features are added or old

problems are repaired. The software is very modular and allows new classes to be added for handling other data formats easily. FabIO and its source-code are freely available to everyone on-line, licensed under the GNU General Public License version 3 (GPLv3). FabIO is also available directly from popular Linux distributions like Debian and Ubuntu.

1.2.2 Implementation

The main language used in the development of FabIO is Python; however, some image formats are compressed and require compression algorithms for reading and writing data. When such algorithms could not be implemented efficiently using Python or NumPy native modules were developed, in i.e. standard C code callable from Python (sometimes generated using Cython). This code has to be compiled for each computer architecture and offers excellent performance. FabIO is only dependent on the NumPy module and has extra features if two other optional Python modules are available. For reading XML files (that are used in EDNA) the Lxml module is required and the Python Image Library, PIL is needed for producing a PIL image for displaying the image in graphical user interfaces and several image-processing operations that are not re-implemented in FabIO. A variety of useful image processing is also available in the `scipy.ndimage` module and in `scikits-image`.

Images can also be displayed in a convenient interactive manner using matplotlib and an IPython shell, which is mainly used for developing data analysis algorithms. Reading and writing procedure of the various TIFF formats is based on the TiffIO code from PyMCA.

In the Python shell, the *fabio* module must be imported prior to reading an image in one of the supported file formats (see Table *Supported formats*, hereafter). The *fabio.open* function creates an instance of the Python class *fabioimage*, from the name of a file. This instance, named *img* hereafter, stores the image data in *img.data* as a 2D NumPy array. Often the image file contains more information than just the intensities of the pixels, e.g. information about how the image is stored and the instrument parameters at the time of the image acquisition, these metadata are usually stored in the file header. Header information, are available in *img.header* as a Python dictionary where keys are strings and values are usually strings or numeric values.

Information in the header about the binary part of the image (compression, endianness, shape) are interpreted however, other metadata are exposed as they are recorded in the file. FabIO allows the user to modify and, where possible, to save this information (the table *Supported formats* summarizes writable formats). Automatic translation between file-formats, even if desirable, is sometimes impossible because not all format have the capability to be extended with additional metadata. Nevertheless FabIO is capable of converting one image data-format into another by taking care of the numerical specifics: for example float arrays are converted to integer arrays if the output format only accepts integers.

1.2.3 FabIO methods

One strength of the implementation in an object oriented language is the possibility to combine functions (or methods) together with data appropriate for specific formats. In addition to the header information and image data, every *fabioimage* instance (returned by *fabio.open*) has methods inherited from *fabioimage* which provide information about the image minimum, maximum and mean values. In addition there are methods which return the file number, name etc. Some of the most important methods are specific for certain formats because the methods are related to how frames in a sequence are handled; these methods are *img.next()*, *img.previous()*, and *img.getframe(n)*. The behaviour of such methods varies depending on the image format: for single-frame format (like `mar345`), *img.next()* will return the image in next file; for multi-frame format (like `GE`), *img.next()* will return the next frame within the same file. For formats which are possibly multi-framed like `EDF`, the behaviour depends on the actual number of frames per file (accessible via the *img.nframes* attribute).

1.3 Usage

1.3.1 Examples

In this section we have collected some basic examples of how FabIO can be employed.

Opening an image:

```
import fabio
im100 = fabio.open('Quartz_0100.tif') # Open image file
print(im0.data[1024,1024])           # Check a pixel value
im101 = im100.next()                 # Open next image
im270 = im1.getframe(270)            # Jump to file number 270: Quartz_0270.tif
```

Normalising the intensity to a value in the header:

```
img = fabio.open('exampleimage0001.edf')
print(img.header)
{'ByteOrder': 'LowByteFirst',
 'DATE (scan begin)': 'Mon Jun 28 21:22:16 2010',
 'ESRFCurrent': '198.099',
 ...
}
# Normalise to beam current and save data
srcur = float(img.header['ESRFCurrent'])
img.data *= 200.0/srcur
img.write('normed_0001.edf')
```

Interactive viewing with matplotlib:

```
from matplotlib import pyplot        # Load matplotlib
pyplot.imshow(img.data)              # Display as an image
pyplot.show()                        # Show GUI window
```

1.4 Future and perspectives

The Hierarchical Data Format version 5 (*hdf5*) is a data format which is increasingly popular for storage of X-ray and neutron data. To name a few facilities the synchrotron Soleil and the neutron sources ISIS, SNS and SINQ already use HDF extensively through the NeXus format. For now, mainly processed or curated data are stored in this format but new detectors are rumoured to provide native output in HDF5. FabIO will rely on H5Py, which already provides a good HDF5 binding for Python, as an external dependency, to be able to read and write such HDF5 files.

In the near future FabIO will be upgraded to work with Python3 (a new version of Python); this change of version will affect some internals FabIO as string and file handling have been altered. This change is already ongoing as many parts of native code in C have already been translated into Cython to smoothe the transition, since Cython generates code compatible with Python3. This also makes it easier to retain backwards compatibility with the earlier Python versions.

1.5 Conclusion

FabIO gives an easy way to read and write 2D images when using the Python computer language. It was originally developed for X-ray diffraction data but now gives an easy way for scientists to access and manipulate their data from a wide range of 2D X-ray detectors. We welcome contributions to further improve the code and hope to add more file formats in the future as well as port the existing code base to the emerging Python3.

1.5.1 Acknowledgements

We acknowledge Andy Götz and Kenneth Evans for extensive testing when including the FabIO reader in the Fable image viewer (Götz et al., 2007). We also thank V. Armando Solé for assistance with his TiffIO reader and Carsten Gundlach for deployment of FabIO at the beamlines i711 and i811, MAX IV, and providing bug reports. We finally acknowledge our colleagues who have reported bugs and helped to improve FabIO. Financial support was granted by the EU 6th Framework NEST/ADVENTURE project TotalCryst (Poulsen et al., 2006).

1.5.2 Citation

Knudsen, E. B., Sørensen, H. O., Wright, J. P., Goret, G. & Kieffer, J. (2013). J. Appl. Cryst. 46, 537-539.

<http://dx.doi.org/10.1107/S0021889813000150>

1.5.3 List of file formats that FabIO can read and write

In alphabetical order. The listed filename extensions are typical examples. FabIO tries to deduce the actual format from the file itself and only uses extensions as a fallback if that fails.

Table 1.1: Supported formats

Python Module	Detector / Format	Extension	Read	Multi-image	Write
ADSC	ADSC Quantum	.img	Yes	No	Yes
Bruker	Bruker formats	.sfrm	Yes	No	Yes
DM3	Gatan Digital Micrograph	.dm3	Yes	No	No
EDF	ESRF data format	.edf	Yes	Yes	Yes
EDNA-XML	Used by EDNA	.xml	Yes	No	No
CBF	CIF binary files	.cbf	Yes	No	Yes
kcd	Nonius KappaCCD	.kccd	Yes	No	No
fit2d mask	Used by Fit2D	.msk	Yes	No	Yes
fit2d spreadsheet	Used by Fit2D	.spr	Yes	No	Yes
GE	General Electric	No	Yes	Yes	No
HiPiC	Hamamatsu CCD	.tif	Yes	No	No
HDF5	Hierarchical data format	.h5	Yes	No	No
marccd	MarCCD/Mar165	.mccd	Yes	No	No
mar345	Mar345 image plate	.mar3450	Yes	No	Yes
OXD	Oxford Diffraction	.img	Yes	No	Yes
Pixi	pixi	.	Yes	No	No
pilatus	Dectris Pilatus Tiff	.tif	Yes	No	Yes
PNM	Portable aNy Map	.pnm	Yes	No	No
Raxis	Rigaku Saxes format	.img	Yes	No	No
TIFF	Tagged Image File Format	.tif	Yes	No	Yes

1.5.4 Adding new file formats

We hope it will be relatively easy to add new file formats to fabio in the future. The basic idea is the following:

1. inherit from `fabioimage` overriding the methods `_readheader`, `read` and optionally `write`. Name your new module `XXXimage` where `XXX` means something (eg `tifimage`).
2. `readheader` fills in a dictionary of “name”:”value” pairs in `self.header`. No one expects to find anything much in there.

3. read fills in self.data with a numpy array holding the image. Some redundant info which also appears are self.dim1 and self.dim2: the image dimensions, self.bpp is the bytes per pixel and self.bytecode is the numpy.dtype.type of the data.
4. The member variables “_need_a_seek_to_read” and “_need_a_real_file” are there in case you have trouble with the transparent handling of bz2 and gz files.
5. Register the file type (extension naming) in `fabioutils.py:FILETYPES`
6. Add your new module as an import into `fabio.openimage`
7. Fill out the magic numbers for your format in `fabio.openimage` if you know them (the characteristic first few bytes in the file)
8. Upload a testimage to the file release system and create a unittest testcase which opens an example of your new format, confirming the image has actually been read in successfully (eg check the mean, max, min and esd are all correct, perhaps orientation too)
9. Run pylint on your code and then please go clean it up. Have a go at mine while you are at it.
10. Bask in the warm glow of appreciation when someone unexpectedly learns they don't need to convert their data into another format

INSTALLATION

FabIO can, as any Python module, be installed from its sources, available on sourceforge but we advice to use binary packages provided for the most common platforms on sourceforge: Windows, MacOSX and Linux. Moreover FabIO is part of the common Linux distributions Ubuntu (since 11.10) and Debian7 where the package is named python-fabio and can be installed via:

```
# apt-get install python-fabio
```

If you are using MS Windows or MacOSX; binary version have been packaged. Windows installers are executable, just download the one corresponding to you python version and run it. MacOSX builds are zipped: unzip them at the right place.

2.1 Dependencies

- Python 2.6 or later (python 3.x is not yet ready)
- numpy - <http://www.numpy.org>

For full functionality of Fabio the following modules need to be installed:

- PIL (python imaging library) - <http://www.pythonware.com>
- lxml (library for reading XSDImages)

2.2 Installation from sources

FabIO can be downloaded from the fable download page on sourceforge.net. Presently the source code has been distributed as a zip package and a compressed tarball. Download either one and unpack it.

```
http://sourceforge.net/projects/fable/files/fabio/
```

e.g.

```
tar xvzf fabio-0.1.4.tar.gz
```

or

```
unzip fabio-0.1.4.zip
```

all files are unpacked into the directory fabio-0.1.4. To install these do

```
cd fabio-0.1.4
```

and install fabio with

```
python setup.py build
sudo python setup.py install
```

most likely you will need to gain root privileges (with sudo in front of the command) to install the built package.

2.3 Development versions

The newest development version can be obtained by checking it out from the subversion (SVN) repository:

```
svn checkout https://svn.sourceforge.net/svnroot/fable/fabio/trunk fabio
cd fabio
python setup.py build
sudo python setup.py install
```

For Ubuntu/Debian users, you will need:

- python-imaging
- python-imaging-tk
- python-numpy
- python-dev

```
sudo apt-get install python-imaging python-imaging-tk python-numpy
```

We provide also a debian-package builder based on stdeb:

```
sudo apt-get install python-stdeb
./build-deb.sh
```

which builds a debian package and installs it in a single command. Handy for testing.

2.4 Test suite

FabIO has a comprehensive test-suite to ensure non regression. When you run the test for the first time, many test images will be download and converted into various compressed format like gzip and bzip2 (this takes a lot of time).

Be sure you have an internet connection and your environment variable `http_proxy` is correctly set-up. For example if you are behind a firewall/proxy:

```
export http_proxy=http://proxy.site.org:3128
```

To run the test:

```
python setup.py test
```

Many tests are there to deal with malformed files, don't worry if the programs complains in warnings about "bad files", it is done on purpose to ensure robustness in FabIO.

FabIO comes with 25 test-suites (110 tests in total) representing a coverage of 60%. This ensures both non regression over time and ease the distribution under different platforms: FabIO runs under Linux, MacOSX and Windows (in each case in 32 and 64 bits) with python versions 2.6 and 2.7.

Table 2.1: Test suite coverage

Name	Stmts	Exec	Cover
fabio/GEimage	94	48	51%
fabio/HiPiCimage	55	7	12%
fabio/OXDimage	285	271	95%
fabio/TiffIO	794	534	67%
fabio/__init__	15	15	100%
fabio/adscimage	79	37	46%
fabio/binaryimage	50	15	30%
fabio/bruker100image	60	13	21%
fabio/brukerimage	212	171	80%
fabio/cbfimage	441	219	49%
fabio/compression	223	136	60%
fabio/converters	17	14	82%
fabio/dm3image	133	16	12%
fabio/edfimage	596	397	66%
fabio/fabioimage	306	193	63%
fabio/fabioutils	322	256	79%
Continued on next page			

Table 2.1 – continued from previous page

fabio/file_series	140	61	43%
fabio/fit2dmaskimage	75	71	94%
fabio/fit2ds spreadsheetimage	47	7	14%
fabio/hdf5image	146	25	17%
fabio/kcdimage	80	65	81%
fabio/mar345image	244	215	88%
fabio/marccdimage	63	56	88%
fabio/mrcimage	96	0	0%
fabio/openimage	104	69	66%
fabio/pilatusimage	34	5	14%
fabio/pixiimage	95	22	23%
fabio/pnmimage	109	21	19%
fabio/raxisimage	98	88	89%
fabio/readbytestream	26	5	19%
fabio/tifimage	167	60	35%
fabio/xsdimage	94	68	72%

CHANGELOG

3.1 From FabIO-0.1.3 to FabIO-0.1.4:

- Work on compatibility with Python3
- Specific debian support with test images included but no auto-generated files
- Image viewer (fabio_viewer) based on Qt4 (Thanks for Gaël Goret)
- Reading images from HDF5 datasets
- Read support for “MRC” images
- Read support for “Pixi detector (Thanks for Jon Wright)
- Read support for “Raxis” images from Rigaku (Thanks to Brian Pauw)
- Write support for fit2d mask images
- Drop support for python 2.5 + Cythonization of other algorithms

3.2 From FabIO-0.1.2 to FabIO-0.1.3:

- Fixed a memory-leak in mar345 module
- Improved support for bruker format (writer & reader)
- Fixed a bug in EDF headers (very long headers)
- Provide template for new file-formats
- Fix a bug related to PIL in new MacOSX
- Allow binary-images to be read from end

3.3 From FabIO-0.1.1 to FabIO-0.1.2:

- Fixed a bug in fabioimage.write (impacted all writers)
- added Sphinx documentation “python setup.py build_doc”
- PyLint compliance of some classes (rename, ...)
- tests from installer with “python setup.py build test”

3.4 From FabIO-0.1.0 to FabIO-0.1.1:

- Merged Mar345 image reader and writer with cython bindings (towards python3 compliance)
- Improve CBF image writing under windows
- Bz2, Gzip and Flat files are managed through a common way ... classes are more (python v2.5) or less (python v2.7) overloaded
- Fast EDF reading if one assumes offsets are the same between files, same for ROIs

3.5 From FabIO-0.0.8 to FabIO-0.1.0:

- OXD reader improved and writer implemented
- Mar345 reader improved and writer implemented
- CBF writer implemented
- Clean-up of the code & bug fixes
- Move towards python3
- Make PIL optional dependency

Python3 is not yet tested but some blocking points have been identified and some fixed.

3.6 From FabIO-0.0.7 to FabIO-0.0.8:

- Support for Tiff using TiffIO module from V.A.Solé
- Clean-up of the code & bug fixes

3.7 From FabIO-0.0.6 to FabIO-0.0.7:

- Support for multi-frames EDF files
- Support for XML images/2D arrays used in EDNA
- new method: `fabio.open(filename)` that is an alias for `fabio.openimage.openimage(filename)`

3.8 From FabIO-0.0.4 to FabIO-0.0.6:

- Support for CBF files from Pilatus detectors
- Support for KCD files from Nonius Kappa CCD images
- write EDF with their native data type (instead of uint16 by default)

FABIO PACKAGE

4.1 fabio Package

FabIO module

4.2 fabio.fabioimage Module

Authors: Henning O. Sorensen & Erik Knudsen Center for Fundamental Research: Metal Structures in Four Dimensions Risoe National Laboratory Frederiksborgvej 399 DK-4000 Roskilde email:erik.knudsen@risoe.dk and Jon Wright, Jerome Kieffer: ESRF

class `fabio.fabioimage.fabioimage` (*data=None, header=None*)

Bases: `object`

A common object for images in fable Contains a numpy array (.data) and dict of meta data (.header)

add (*other*)

Add another Image - warning, does not clip to 16 bit images by default

static checkData (*data=None*)

Empty for fabioimage but may be populated by others classes, especially for format accepting only integers

static checkHeader (*header=None*)

Empty for fabioimage but may be populated by others classes

classname

Retrieves the name of the class :return: the name of the class

convert (*dest*)

Convert a fabioimage object into another fabioimage object (with possible conversions) :param dest: destination type "EDF", "edfimage" or the class itself

getclassname ()

Retrieves the name of the class :return: the name of the class

getframe (*num*)

returns the file numbered 'num' in the series as a fabioimage

getheader ()

returns self.header

getmax ()

Find max value in self.data, caching for the future

getmean ()
 return the mean

getmin ()
 Find min value in self.data, caching for the future

getstddev ()
 return the standard deviation

integrate_area (*coords*)
 Sums up a region of interest if len(coords) == 4 -> convert coords to slices if len(coords) == 2 -> use as slices floor -> ? removed as unused in the function.

load (**arg*, ***kwarg*)
 Wrapper for read

make_slice (*coords*)
 Convert a len(4) set of coords into a len(2) tuple (pair) of slice objects the latter are immutable, meaning the roi can be cached

next ()
 returns the next file in the series as a fabioimage

previous ()
 returns the previous file in the series as a fabioimage

read (*filename*, *frame=None*)
 To be overridden - fill in self.header and self.data

readROI (*filename*, *frame=None*, *coords=None*)
 Method reading Region of Interest. This implementation is the trivial one, just doing read and crop

readheader (*filename*)
 Call the _readheader function...

rebin (*x_rebin_fact*, *y_rebin_fact*, *keep_I=True*)
 Rebin the data and adjust dims :param x_rebin_fact: x binning factor :param y_rebin_fact: y binning factor :param keep_I: shall the signal increase ? :type x_rebin_fact: int :type y_rebin_fact: int :type keep_I: boolean

resetvals ()
 Reset cache - call on changing data

save (*fname*)
 wrapper for write

toPIL16 (*filename=None*)
 Convert to Python Imaging Library 16 bit greyscale image
 FIXME - this should be handled by the libraries now

update_header (***kws*)
 update the header entries by default pass in a dict of key, values.

write (*fname*)
 To be overwritten - write the file

fabio.fabioimage.test ()
 check some basic fabioimage functionality

4.3 fabio.fabioutils Module

General purpose utilities functions for fabio

class fabio.fabioutils.**BZ2File** (*name, mode='r', buffering=0, compresslevel=9*)
 Bases: bz2.BZ2File

Wrapper with lock

getSize()

setSize (*value*)

size

class fabio.fabioutils.**DebugSemaphore** (**arg, **kwarg*)
 Bases: threading._Semaphore

threading.Semaphore like class with helper for fighting dead-locks

acquire (**arg, **kwarg*)

release (**arg, **kwarg*)

class fabio.fabioutils.**File** (*name, mode='rb', buffering=0*)
 Bases: file

wrapper for “file” with locking

getSize()

setSize (*size*)

size

class fabio.fabioutils.**FilenameObject** (*stem=None, num=None, directory=None, format=None, extension=None, postnum=None, digits=4, filename=None*)

Bases: object

The ‘meaning’ of a filename ...

deconstruct_filename (*filename*)

Break up a filename to get image type and number

str()

Return a string representation

tostring()

convert yourself to a string

class fabio.fabioutils.**GzipFile** (*filename=None, mode=None, compresslevel=9, fileobj=None*)
 Bases: gzip.GzipFile

Just a wrapper for gzip.GzipFile providing the correct seek capabilities for python 2.5

closed

getSize()

seek (*offset, whence=0*)

Move to new file position.

Argument offset is a byte count. Optional argument whence defaults to 0 (offset from start of file, offset should be >= 0); other values are 1 (move relative to current position, positive or negative), and 2 (move relative to end of file, usually negative, although many platforms allow seeking beyond the end of a file).

If the file is opened in text mode, only offsets returned by `tell()` are legal. Use of other offsets causes undefined behavior.

This is a wrapper for `seek` to ensure compatibility with old python 2.5

Warning: `Seek from end` is not supported (works only for single blocks !!!) This implements a hack

setSize (*value*)

size

class `fabio.fabioutils.StringIO` (*data, fname=None, mode='r'*)

Bases: `StringIO.StringIO`

just an interface providing the name and mode property to a `StringIO`

BugFix for MacOSX mainly

getSize ()

setSize (*size*)

size

class `fabio.fabioutils.UnknownCompressedFile` (*name, mode='rb', buffering=0*)

Bases: `fabio.fabioutils.File`

wrapper for “File” with locking

`fabio.fabioutils.construct_filename` (*filename, frame=None*)

Try to construct the filename for a given frame

`fabio.fabioutils.deconstruct_filename` (*filename*)

Function for backward compatibility. Deprecated

`fabio.fabioutils.deprecated` (*func*)

used to deprecate a function/method: prints a lot of warning messages to enforce the modification of the code

`fabio.fabioutils.extract_filenumber` (*name*)

extract file number

`fabio.fabioutils.getnum` (*name*)

try to figure out a file number # guess it starts at the back

`fabio.fabioutils.isAscii` (*name, listExcluded=None*)

Parameters

- **name** – string to check
- **listExcluded** – list of char or string excluded.

Returns True or False whether name is pure ascii or not

`fabio.fabioutils.jump_filename` (*name, num, padding=True*)

jump to number

`fabio.fabioutils.next_filename` (*name, padding=True*)

increment number

`fabio.fabioutils.nice_int` (*s*)

Workaround that `int('1.0')` raises an exception

Parameters *s* – string to be converted to integer

`fabio.fabioutils.numstem` (*name*)

can't see how to do without reversing strings Match 1 or more digits going backwards from the end of the string

```
fabio.fabioutils.pad(mystr, pattern=' ', size=80)
    Performs the padding of the string to the right size with the right pattern

fabio.fabioutils.previous_filename(name, padding=True)
    decrement number

fabio.fabioutils.toAscii(name, excluded=None)
```

Parameters

- **name** – string to check
- **excluded** – tuple of char or string excluded (not list: they are mutable).

Returns the name with all non valid char removed

4.4 fabio.file_series Module

4.4.1 Authors:

- Henning O. Sorensen & Erik Knudsen Center for Fundamental Research: Metal Structures in Four Dimensions
Risoe National Laboratory Frederiksborgvej 399 DK-4000 Roskilde email:erik.knudsen@risoe.dk
- Jon Wright, ESRF

```
class fabio.file_series.file_series(list_of_strings)
    Bases: list
```

Represents a series of files to iterate has an idea of a current position to do next and prev

You also get from the list python superclass: append count extend insert pop remove reverse sort

```
current()
    Current position in a sequence
```

```
current_image()
    Current image in sequence
```

Returns fabioimage

```
current_object()
    Current image in sequence
```

Returns file_object

```
first()
    First image in series
```

```
first_image()
    First image in a sequence
```

Returns fabioimage

```
first_object()
    First image in a sequence
```

Returns file_object

```
jump(num)
    Goto a position in sequence
```

```
jump_image(num)
    Jump to and read image
```

Returns fabioimage

jump_object (*num*)
 Jump to and read image

Returns file_object

last ()
 Last in series

last_image ()
 Last image in a sequence

Returns fabioimage

last_object ()
 Last image in a sequence

Returns file_object

len ()
 Number of files

next ()
 Next in a sequence

next_image ()
 Return the next image

Returns fabioimage

next_object ()
 Return the next image

Returns file_object

previous ()
 Prev in a sequence

previous_image ()
 Return the previous image

Returns fabioimage

previous_object ()
 Return the previous image

Returns file_object

class fabio.file_series.**filename_series** (*filename*)
 Much like the others, but created from a string filename

current ()
 return current filename string

current_image ()
 returns the current image as a fabioimage

current_object ()
 returns the current filename as a fabio.FileNameObject

jump (*num*)
 jump to a specific number

jump_image (*num*)
 returns the image number as a fabioimage

jump_object (*num*)
returns the filename num as a fabio.FilenameObject

next ()
increment number

next_image ()
returns the next image as a fabioimage

next_object ()
returns the next filename as a fabio.FilenameObject

prev_image ()
returns the previos image as a fabioimage

previous ()
decrement number

previous_object ()
returns the previous filename as a fabio.FilenameObject

`fabio.file_series.new_file_series` (*first_object*, *nimages=0*, *step=1*, *traceback=False*)

A generator function that creates a file series starting from a a fabioimage. Iterates through all images in a file (if more than 1), then proceeds to the next file as determined by `fabio.next_filename`.

Parameters

- **first_object** – the starting fabioimage, which will be the first one yielded in the sequence
- **nimages** – the maximum number of images to consider step: step size, will yield the first and every step'th image until nimages is reached. (e.g. nimages = 5, step = 2 will yield 3 images (0, 2, 4))
- **traceback** – if True causes it to print a traceback in the event of an exception (missing image, etc.). Otherwise the calling routine can handle the exception as it chooses
- **yields** – the next fabioimage in the series. In the event there is an exception, it yields the `sys.exec_info` for the exception instead. `sys.exec_info` is a tuple: (`exceptionType`, `exceptionValue`, `exceptionTraceback`) from which all the exception information can be obtained.

Suggested usage:

```
for obj in new_file_series( ... ):
    if not isinstance(obj, fabio.fabioimage.fabioimage ):
        # deal with errors like missing images, non readable files, etc
        # e.g.
        traceback.print_exception(obj[0], obj[1], obj[2])
```

`fabio.file_series.new_file_series0` (*first_object*, *first=None*, *last=None*, *step=1*)

Created from a fabio image first and last are file numbers

class `fabio.file_series.numbered_file_series` (*stem*, *first*, *last*, *extension*, *digits=4*,
padding='Y', *step=1*)

Bases: `fabio.file_series.file_series`

`mydata0001.edf` = “mydata” + 0001 + “.edf” `mydata0002.edf` = “mydata” + 0002 + “.edf” `mydata0003.edf` = “mydata” + 0003 + “.edf”

4.5 fabio.openimage Module

Authors: Henning O. Sorensen & Erik Knudsen Center for Fundamental Research: Metal Structures in Four Dimensions Risoe National Laboratory Frederiksborgvej 399 DK-4000 Roskilde email:henning.sorensen@risoe.dk

mods for fabio by JPW

`fabio.openimage.do_magic` (*byts*)
Try to interpret the bytes starting the file as a magic number

`fabio.openimage.openheader` (*filename*)
return only the header

`fabio.openimage.openimage` (*filename, frame=None*)
Try to open an image

4.6 fabio.adscimage Module

Authors: Henning O. Sorensen & Erik Knudsen Center for Fundamental Research: Metal Structures in Four Dimensions Risoe National Laboratory Frederiksborgvej 399 DK-4000 Roskilde email:erik.knudsen@risoe.dk

- mods for fabio by JPW

class `fabio.adscimage.adscimage` (**args, **kwargs*)
Bases: `fabio.fabioimage.fabioimage`

Read an image in ADSC format (quite similar to edf?)

read (*fname, frame=None*)
read in the file

write (*fname*)
Write adsc format

`fabio.adscimage.test` ()
testcase

4.7 fabio.binaryimage Module

Authors: Gael Goret, Jerome Kieffer, ESRF, France Emails: gael.goret@esrf.fr, jerome.kieffer@esrf.fr

Brian Richard Pauw <brian@stack.nl>

Binary files images are simple none-compressed 2D images only defined by their : data-type, dimensions, byte order and offset

This simple library has been made for manipulating exotic/unknown files format.

class `fabio.binaryimage.binaryimage` (**args, **kwargs*)
Bases: `fabio.fabioimage.fabioimage`

This simple library has been made for manipulating exotic/unknown files format.

Binary files images are simple none-compressed 2D images only defined by their : data-type, dimensions, byte order and offset

if offset is set to a negative value, the image is read using the last data but n data in the file, skipping any header.

estimate_offset_value (*fname, dim1, dim2, bytecode='int32'*)

Estimates the size of a file

read (*fname, dim1, dim2, offset=0, bytecode='int32', endian='<'*)

Read a binary image

Parameters

- **fname** (*string*) – file name
- **dim1** – image dimensions (Fast index)
- **dim2** – image dimensions (Slow index)
- **offset** – starting position of the data-block. If negative, starts at the end.
- **bytecode** – can be “int8”, “int16”, “int32”, “int64”, “uint8”, “uint16”, “uint32”, “uint64”, “float32”, “float64”, ...
- **endian** – among short or long endian (“<” or “>”)

static swap_needed (*endian*)

Decide if we need to byteswap

write (*fname*)

4.8 fabio.bruker100image Module

class fabio.bruker100image.**bruker100image** (*data=None, header=None*)

Bases: fabio.brukerimage.brukerimage

read (*fname, frame=None*)

toPIL16 (*filename=None*)

4.9 fabio.brukerimage Module

Authors: Henning O. Sorensen & Erik Knudsen Center for Fundamental Research: Metal Structures in Four Dimensions Risoe National Laboratory Frederiksborgvej 399 DK-4000 Roskilde email:erik.knudsen@risoe.dk

Based on: **openbruker, readbruker, readbrukerheader functions in the opendata** module of ImageD11 written by Jon Wright, ESRF, Grenoble, France

Writer by Jérôme Kieffer, ESRF, Grenoble, France

class fabio.brukerimage.**brukerimage** (*data=None, header=None*)

Bases: fabio.fabioimage.fabioimage

Read and eventually write ID11 bruker (eg smart6500) images

TODO: int32 -> float32 conversion according to the “linear” keyword. This is done and works but we need to check with other program that we are applying the right formula and not the reciprocal one.

basic_translate (*fname=None*)

Does some basic population of the headers so that the writing is possible

calc_bpp (*data=None, max_entry=4096*)

Calculate the number of byte per pixel to get an optimal overflow table.

Returns byte per pixel

```

gen_header ()
    Generate headers (with some magic and guesses) :param format can be 86 or 100

gen_overflow ()
    Generate an overflow table

read (fname, frame=None)
    Read in and unpack the pixels (including overflow table)

write (fname)
    Write a bruker image

fabio.brukerimage.test ()
    a testcase

```

4.10 fabio.cbfimage Module

Authors: Jérôme Kieffer, ESRF email:jerome.kieffer@esrf.fr

Cif Binary Files images are 2D images written by the Pilatus detector and others. They use a modified (simplified) byte-offset algorithm.

CIF is a library for manipulating Crystallographic information files and tries to conform to the specification of the IUCR

```

class fabio.cbfimage.CIF (_strFilename=None)
    Bases: dict

    This is the CIF class, it represents the CIF dictionary; and as a python dictionary thus inherits from the dict
    built in class.

    static LoopHasKey (loop, key)
        Returns True if the key (string) exist in the array called loop

    exists (sKey)
        Check if the key exists in the CIF and is non empty. :param sKey: CIF key :type sKey: string :param cif:
        CIF dictionary :return: True if the key exists in the CIF dictionary and is non empty :rtype: boolean

    existsInLoop (sKey)
        Check if the key exists in the CIF dictionary. :param sKey: CIF key :type sKey: string :param cif: CIF
        dictionary :return: True if the key exists in the CIF dictionary and is non empty :rtype: boolean

    static isAscii (_strIn)
        Check if all characters in a string are ascii,

        Parameters _strIn (python string) – input string

        Returns boolean

        Return type boolean

    loadCHIPLLOT (_strFilename)
        Load the powder diffraction CHIPLLOT file and returns the pd_CIF dictionary in the object

        Parameters _strFilename (string) – the name of the file to open

        Returns the CIF object corresponding to the powder diffraction

        Return type dictionary

```

loadCIF (*_strFilename*, *_bKeepComment=False*)

Load the CIF file and populates the CIF dictionary into the object :param *_strFilename*: the name of the file to open :type *_strFilename*: string :param *_strFilename*: the name of the file to open :type *_strFilename*: string :return: None

pop (*key*)

popitem (*key*)

readCIF (*_strFilename*, *_bKeepComment=False*)

Load the CIF file and populates the CIF dictionary into the object :param *_strFilename*: the name of the file to open :type *_strFilename*: string :param *_strFilename*: the name of the file to open :type *_strFilename*: string :return: None

saveCIF (*_strFilename='test.cif'*, *linesep='n'*, *binary=False*)

Transforms the CIF object in string then write it into the given file :param *_strFilename*: the of the file to be written :param *linesep*: line separation used (to force compatibility with windows/unix) :param *binary*: Shall we write the data as binary (True only for imageCIF/CBF) :type *param*: string

tostring (*_strFilename=None*, *linesep='n'*)

Converts a cif dictionary to a string according to the CIF syntax

Parameters *_strFilename* (*string*) – the name of the filename to be appended in the header of the CIF file

Returns a sting that corresponds to the content of the CIF - file.

class `fabio.cbimage.cbimage` (*data=None*, *header=None*, *fname=None*)

Bases: `fabio.fabioimage.fabioimage`

Read the Cif Binary File data format

static **checkData** (*data=None*)

read (*fname*, *frame=None*)

Read in header into self.header and the data into self.data

write (*fname*)

write the file in CBF format :param *fname*: name of the file :type: string

4.11 fabio.dm3image Module

Authors: Henning O. Sorensen & Erik Knudsen

Center for Fundamental Research: Metal Structures in Four Dimensions Risoe National Laboratory
Frederiksborgvej 399 DK-4000 Roskilde email:erik.knudsen@risoe.dk

- Jon Wright, ESRF

class `fabio.dm3image.dm3image` (**args*, ***kwargs*)

Bases: `fabio.fabioimage.fabioimage`

Read and try to write the dm3 data format

read (*fname*, *frame=None*)

read_data ()

read_tag_entry ()

read_tag_group ()

```
read_tag_type ()
readbytes (bytes_to_read, format, swap=True)
```

4.12 fabio.edfimage Module

License: GPLv2+

4.12.1 Authors:

- Henning O. Sorensen & Erik Knudsen: Center for Fundamental Research: Metal Structures in Four Dimensions; Risoe National Laboratory; Frederiksborgvej 399; DK-4000 Roskilde; email:erik.knudsen@risoe.dk
- Jon Wright & Jérôme Kieffer: European Synchrotron Radiation Facility; Grenoble (France)

```
class fabio.edfimage.Frame (data=None, header=None, header_keys=None, number=None)
    Bases: object
```

A class representing a single frame in an EDF file

bytecode

data

Unpack a binary blob according to the specification given in the header

Returns dataset as numpy.ndarray

getByteCode ()

getData ()

Unpack a binary blob according to the specification given in the header

Returns dataset as numpy.ndarray

getEdfBlock (*force_type=None, fit2dMode=False*)

Parameters

- **force_type** (*string or numpy.dtype*) – type of the dataset to be enforced like “float64” or “uint16”
- **fit2dMode** (*boolean*) – enforce compatibility with fit2d and starts counting number of images at 1

Returns ascii header block

Return type python string with the concatenation of the ascii header and the binary data block

parseheader (*block*)

Parse the header in some EDF format from an already open file

Parameters **block** (*string, should be full ascii*) – string representing the header block

Returns size of the binary blob

setByteCode (*_iVal*)

setData (*npa=None*)

Setter for data in edf frame

swap_needed ()

Decide if we need to byteswap

```

class fabio.edfimage.edfimage (data=None, header=None, header_keys=None, frames=None)
    Bases: fabio.fabioimage.fabioimage

    Read and try to write the ESRF edf data format

    appendFrame (frame=None, data=None, header=None)
        Method used add a frame to an EDF file :param frame: frame to append to edf image :type frame: instance
        of Frame :return: None

    bpp

    bytecode

    capsHeader
        property: capsHeader of EDF file, i.e. the keys of the header in UPPER case.

    static checkHeader (header=None)
        Empty for fabioimage but may be populated by others classes

    data
        property: data of EDF file

    delCapsHeader ()
        deleter for edf capsHeader

    delData ()
        deleter for edf Data

    delHeader ()
        Deleter for edf header

    delHeaderKeys ()
        Deleter for edf header_keys

    deleteFrame (frameNb=None)
        Method used to remove a frame from an EDF image. by default the last one is removed. :param frameNb:
        frame number to remove, by default the last. :type frameNb: integer :return: None

    dim1

    dim2

    dims

    fastReadData (filename=None)
        This is a special method that will read and return the data from another file ... The aim is performances, ...
        but only supports uncompressed files.

        Returns data from another file using positions from current edfimage

    fastReadROI (filename, coords=None)
        Method reading Region of Interest of another file based on metadata available in current edfimage. The
        aim is performances, ... but only supports uncompressed files.

        Returns ROI-data from another file using positions from current edfimage

        Return type numpy 2darray

    getBpp ()

    getByteCode ()

    getCapsHeader ()
        getter for edf headers keys in upper case :return: data for current frame :rtype: dict

```

getData ()
 getter for edf Data :return: data for current frame :rtype: numpy.ndarray

getDim1 ()

getDim2 ()

getDims ()

getHeader ()
 Getter for the headers. used by the property header,

getHeaderKeys ()
 Getter for edf header_keys

getNbFrames ()
 Getter for number of frames

getframe (num)
 returns the file numbered 'num' in the series as a fabioimage

header
 property: header of EDF file

header_keys
 property: header_keys of EDF file

next ()
 returns the next file in the series as a fabioimage

nframes
 Getter for number of frames

previous ()
 returns the previous file in the series as a fabioimage

read (fname, frame=None)
 Read in header into self.header and the data into self.data

setBpp (_iVal)

setByteCode (_iVal)

setCapsHeader (_data)
 Enforces the propagation of the header_keys to the list of frames :param _data: numpy array representing data

setData (_data)
 Enforces the propagation of the data to the list of frames :param _data: numpy array representing data

setDim1 (_iVal)

setDim2 (_iVal)

setHeader (_dictHeader)
 Enforces the propagation of the header to the list of frames

setHeaderKeys (_listtHeader)
 Enforces the propagation of the header_keys to the list of frames :param _listtHeader: list of the (ordered) keys in the header :type _listtHeader: python list

setNbFrames (val)
 Setter for number of frames ... should do nothing. Here just to avoid bugs

swap_needed()

Decide if we need to byteswap

unpack()

Unpack a binary blob according to the specification given in the header and return the dataset

Returns dataset as numpy.ndarray

write(*fname*, *force_type=None*, *fit2dMode=False*)

Try to write a file check we can write zipped also mimics that fabian was writing uint16 (we sometimes want floats)

Parameters *force_type* – can be numpy.uint16 or simply “float”

Returns None

4.13 fabio.fit2dmaskimage Module

Author: Andy Hammersley, ESRF Translation into python/fabio: Jon Wright, ESRF. Writer: Jérôme Kieffer

class fabio.fit2dmaskimage.**fit2dmaskimage**(*data=None*, *header=None*)

Bases: fabio.fabioimage.fabioimage

Read and try to write Andy Hammersley’s mask format

static checkData(*data=None*)

read(*fname*, *frame=None*)

Read in header into self.header and the data into self.data

write(*fname*)

Try to write a file

4.14 fabio.fit2dspreadsheetimage Module

Read the fit2d ascii image output

- Jon Wright, ESRF

class fabio.fit2dspreadsheetimage.**fit2dspreadsheetimage**(*data=None*, *header=None*)

Bases: fabio.fabioimage.fabioimage

Read a fit2d ascii format

read(*fname*, *frame=None*)

Read in header into self.header and the data into self.data

4.15 fabio.GEimage Module

class fabio.GEimage.**GEimage**(*data=None*, *header=None*)

Bases: fabio.fabioimage.fabioimage

getframe(*num*)

Returns a frame as a new fabioimage object

```

next ()
    Get the next image in a series as a fabio image

previous ()
    Get the previous image in a series as a fabio image

read (fname, frame=None)
    Read in header into self.header and the data into self.data

write (fname, force_type=<type 'numpy.uint16'>)
    Not yet implemented

fabio.GEImage.demo ()

```

4.16 fabio.hdf5mage Module

HDF5 image for FabIO

Authors: Jerome Kieffer email: Jerome.Kieffer@terre-adelie.org

Specifications: input should being the form:

hdf5://filename:path[slice]

```

class fabio.hdf5image.HDF5location (filename=None, h5path=None, slices=None, url=None)
    Bases: object

    Handle URL like:
    hdf5://filename:path[slice]

    parse (url)
        Analyse a string of the form hdf5://filename:path[slice]

        Parameters url – string of form of an hdf5-url

    set_index (idx)
        Set the current frame to idx

    to_url ()
        convert an HDF5 locate into an URL

class fabio.hdf5image.hdf5image (*arg, **kwargs)
    Bases: fabio.fabioimage.fabioimage

    FabIO image class for Images from an HDF file

    getframe (num)
        Returns a frame as a new fabioimage object :param num: frame number

    next ()
        Get the next image in a series as a fabio image

    previous ()
        Get the previous image in a series as a fabio image

    read (fname, frame=None)
        try to read image :param fname: name of the file as hdf5://filename:path[slice]

    write (fname, force_type=<type 'numpy.uint16'>)

```

4.17 fabio.HiPiCimage Module

Authors: Henning O. Sorensen & Erik Knudsen

Center for Fundamental Research: Metal Structures in Four Dimensions Risoe National Laboratory
Frederiksborgvej 399 DK-4000 Roskilde email:erik.knudsen@risoe.dk

- Jon Wright, ESRF

Information about the file format from Masakatzu Kobayashi is highly appreciated

```
class fabio.HiPiCimage.HiPiCimage (data=None, header=None)
```

```
    Bases: fabio.fabioimage.fabioimage
```

Read HiPic images e.g. collected with a Hamamatsu CCD camera

```
    read (fname, frame=None)
```

Read in header into self.header and the data into self.data

4.18 fabio.kcdimage Module

Authors: Jerome Kieffer, ESRF email:jerome.kieffer@esrf.fr

kcd images are 2D images written by the old KappaCCD diffractometer built by Nonius in the 1990's Based on the edfimage.py parser.

```
class fabio.kcdimage.kcdimage (data=None, header=None)
```

```
    Bases: fabio.fabioimage.fabioimage
```

Read the Nonius kcd data format

```
    static checkData (data=None)
```

```
    read (fname, frame=None)
```

Read in header into self.header and the data into self.data

4.19 fabio.mar345image Module

4.19.1 Authors:

- Henning O. Sorensen & Erik Knudsen: Center for Fundamental Research: Metal Structures in Four Dimensions; Risoe National Laboratory; Frederiksborgvej 399; DK-4000 Roskilde; email:erik.knudsen@risoe.dk
- Jon Wright, Jérôme Kieffer & Gaël Goret: European Synchrotron Radiation Facility; Grenoble (France)

```
class fabio.mar345image.mar345image (*args, **kwargs)
```

```
    Bases: fabio.fabioimage.fabioimage
```

```
    static checkData (data=None)
```

```
    nb_overflow_pixels ()
```

```
    read (fname, frame=None)
```

Read a mar345 image

```
    write (fname)
```

Try to write mar345 file. This is still in beta version. It uses CCP4 (LGPL) PCK1 algo from JPA

4.20 fabio.mrcimage Module

MRC image for FabIO

Authors: Jerome Kieffer email: Jerome.Kieffer@terre-adelie.org

Specifications from: http://ami.scripps.edu/software/mrctools/mrc_specification.php

```
class fabio.mrcimage.mrcimage (*arg, **kwargs)
    Bases: fabio.fabioimage.fabioimage

    FabIO image class for Images from a mrc image stack

    getframe (num)
        Returns a frame as a new fabioimage object :param num: frame number

    next ()
        Get the next image in a series as a fabio image

    previous ()
        Get the previous image in a series as a fabio image

    read (fname, frame=None)
        try to read image :param fname: name of the file :param frame:

    write (fname, force_type=<type 'numpy.uint16'>)
```

4.21 fabio.marccdimage Module

4.21.1 Authors:

- Henning O. Sorensen & Erik Knudsen: Center for Fundamental Research: Metal Structures in Four Dimensions; Risoe National Laboratory; Frederiksborgvej 399; DK-4000 Roskilde; email:erik.knudsen@risoe.dk
- Jon Wright: European Synchrotron Radiation Facility; Grenoble (France)

marccdimage can read MarCCD and MarMosaic images including header info.

JPW : Use a parser in case of typos (sorry?)

```
fabio.marccdimage.interpret_header (header, fmt, names)
    given a format and header interpret it
```

```
fabio.marccdimage.make_format (c_def_string)
    Reads the header definition in c and makes the format string to pass to struct.unpack
```

```
class fabio.marccdimage.marccdimage (*args, **kws)
    Bases: fabio.tifimage.tifimage

    Read in data in mar ccd format, also MarMosaic images, including header info
```

4.22 fabio.OXDImage Module

Reads Oxford Diffraction Sapphire 3 images

4.22.1 Authors:

- Henning O. Sorensen & Erik Knudsen: Center for Fundamental Research: Metal Structures in Four Dimensions; Risoe National Laboratory; Frederiksborgvej 399; DK-4000 Roskilde; email:erik.knudsen@risoe.dk
- Jon Wright, Jérôme Kieffer & Gaël Goret: European Synchrotron Radiation Facility; Grenoble (France)

```
class fabio.OXDImage.OXDImage (data=None, header=None)
    Bases: fabio.fabioimage.fabioimage
    Oxford Diffraction Sapphire 3 images reader/writer class

    static checkData (data=None)

    getCompressionRatio ()
        calculate the compression factor obtained vs raw data

    read (fname, frame=None)

        Read in header into self.header and the data into self.data

    write (fname)
        Write Oxford diffraction images: this is still beta :param fname: output filename

class fabio.OXDImage.Section (size, dictHeader)
    Bases: object
    Small helper class for writing binary headers

    getSize (dtype)

    setData (key, offset, dtype, default=None)

        Parameters
        • offset – int, starting position in the section
        • key – name of the header key
        • dtype – type of the data to insert (defines the size!)
```

4.23 fabio.pilatusimage Module

4.23.1 Authors:

- Henning O. Sorensen & Erik Knudsen: Center for Fundamental Research: Metal Structures in Four Dimensions; Risoe National Laboratory; Frederiksborgvej 399; DK-4000 Roskilde; email:erik.knudsen@risoe.dk
- Jon Wright: European Synchrotron Radiation Facility; Grenoble (France)

```
class fabio.pilatusimage.pilatusimage (*args, **kwargs)
    Bases: fabio.tifimage.tifimage
    Read in Pilatus format, also pilatus images, including header info
```

4.24 fabio.pixiimage Module

Author: Jon Wright, ESRF.

```
fabio.pixiimage.demo (fname)
```

```
class fabio.pixiimage.pixiimage (data=None, header=None)
    Bases: fabio.fabioimage.fabioimage

    getframe (num)
        Returns a frame as a new fabioimage object

    next ()
        Get the next image in a series as a fabio image

    previous ()
        Get the previous image in a series as a fabio image

    read (fname, frame=None)

    write (fname, force_type=<type 'numpy.uint16'>)
```

4.25 fabio.pnmimage Module

Authors: Henning O. Sorensen & Erik Knudsen Center for Fundamental Research: Metal Structures in Four Dimensions Risoe National Laboratory Frederiksborgvej 399 DK-4000 Roskilde email:henning.sorensen@risoe.dk

```
class fabio.pnmimage.pnmimage (*arg, **kwargs)
    Bases: fabio.fabioimage.fabioimage

    P1dec (buf, bytecode)

    P2dec (buf, bytecode)

    static P3dec (buf, bytecode)

    static P4dec (buf, bytecode)

    P5dec (buf, bytecode)

    static P6dec (buf, bytecode)

    static P7dec (buf, bytecode)

    static checkData (data=None)

    read (fname, frame=None)
        try to read PNM images :param fname: name of the file :param frame: not relevant here! PNM is always
        single framed

    write (filename)
```

4.26 fabio.raxisimage Module

Authors: Brian R. Pauw email: brian@stack.nl

Written using information gleaned from the ReadRAXISImage program written by T. L. Hendrixson, made available by Rigaku Americas. Available at: <http://www.rigaku.com/downloads/software/readimage.html>

```
class fabio.raxisimage.raxisimage (*arg, **kwargs)
    Bases: fabio.fabioimage.fabioimage
```

FabIO image class to read Rigaku RAXIS image files. Write functions are not planned as there are plenty of more suitable file formats available for storing detector data. In particular, the MSB used in Rigaku files is used in an uncommon way: it is used as a *multiply-by* flag rather than a normal image value bit. While it is said to

multiply by the value specified in the header, there is at least one case where this is found not to hold, so YMMV and be careful.

```
read (fname, frame=None)
    try to read image :param fname: name of the file :param frame:

rigakuKeys ()

swap_needed ()
    not sure if this function is needed
```

4.27 fabio.tifimage Module

FabIO class for dealing with TIFF images. In facts wraps TiffIO from V. Armando Solé (available in PyMca) or falls back to PIL

4.27.1 Authors:

- Henning O. Sorensen & Erik Knudsen: Center for Fundamental Research: Metal Structures in Four Dimensions; Risoe National Laboratory; Frederiksborgvej 399; DK-4000 Roskilde; email:erik.knudsen@risoe.dk
- Jérôme Kieffer: European Synchrotron Radiation Facility; Grenoble (France)

License: GPLv3+

```
class fabio.tifimage.Image_File_Directory (instring=None, offset=-1)
    Bases: object

    unpack (instring, offset=-1)

class fabio.tifimage.Image_File_Directory_entry (tag=0, tag_type=0, count=0, offset=0)
    Bases: object

    extract_data (full_string)

    unpack (strInput)

class fabio.tifimage.Tiff_header (string)
    Bases: object

class fabio.tifimage.tifimage (*args, **kws)
    Bases: fabio.fabioimage.fabioimage
    Images in TIF format Wraps TiffIO

    read (fname, frame=None)
        Wrapper for TiffIO.

    write (fname)
        Overrides the fabioimage.write method and provides a simple TIFF image writer. :param fname: name of
        the file to save the image to @tag_type fname: string or unicode (file?)...
```

4.28 fabio.xsdimage Module

Authors: Jérôme Kieffer, ESRF email:jerome.kieffer@esrf.fr

XSDImage are XML files containing numpy arrays

class `fabio.xsdimimage.xsdimimage` (*data=None, header=None, fname=None*)
 Bases: `fabio.fabioimage.fabioimage`
 Read the XSDataImage XML File data format
read (*fname, frame=None*)

4.29 fabio.compression Module

Authors: Jérôme Kieffer, ESRF email:jerome.kieffer@esrf.fr

FabIO library containing compression and decompression algorithm for various

`fabio.compression.compByteOffset_numpy` (*data*)
 Compress a dataset into a string using the byte_offset algorithm

Parameters *data* – ndarray

Returns string/bytes with compressed data

test = numpy.array([0,1,2,127,0,1,2,128,0,1,2,32767,0,1,2,32768,0,1,2,2147483647,0,1,2,2147483648,0,1,2,128,129,130,32767,3

`fabio.compression.compPCK` (*data*)
 Modified CCP4 pck compressor used in MAR345 images

Parameters *data* – numpy.ndarray (square array)

Returns compressed stream

`fabio.compression.compTY1` (*data*)
 Modified byte offset compressor used in Oxford Diffraction images

Parameters *data* – numpy.ndarray with the input data (integers!)

Returns 3-tuple of strings: raw_8,raw_16,raw_32 containing raw data with integer of the given size

`fabio.compression.decByteOffset_cython` (*stream, size=None*)

Analyze a stream of char with any length of exception: 2, 4, or 8 bytes integers

Parameters

- **stream** – string representing the compressed data
- **size** – the size of the output array (of longInts)

Returns 1D-ndarray

`fabio.compression.decByteOffset_numpy` (*stream, size=None*)

Analyze a stream of char with any length of exception: 2, 4, or 8 bytes integers

Parameters

- **stream** – string representing the compressed data
- **size** – the size of the output array (of longInts)

Returns 1D-ndarray

`fabio.compression.decByteOffset_python` (*stream, size*)
 Analyze a stream of char with any length of exception (2,4, or 8 bytes integers)

Parameters

- **stream** – string representing the compressed data
- **size** – the size of the output array (of longInts)

Returns 1D-ndarray

`fabio.compression.decByteOffset_weave (stream, size)`

Analyze a stream of char with any length of exception (2,4, or 8 bytes integers)

Parameters

- **stream** – string representing the compressed data
- **size** – the size of the output array (of longInts)

Returns 1D-ndarray

`fabio.compression.decBzip2 (stream)`

Decompress a chunk of data using the bzip2 algorithm from Python

`fabio.compression.decGzip (stream)`

Decompress a chunk of data using the gzip algorithm from Python or alternatives if possible

`fabio.compression.decKM4CCD (raw_8, raw_16=None, raw_32=None)`

Modified byte offset decompressor used in Oxford Diffraction images

Parameters

- **raw_8** – strings containing raw data with integer 8 bits
- **raw_16** – strings containing raw data with integer 16 bits
- **raw_32** – strings containing raw data with integer 32 bits

Returns numpy.ndarray

`fabio.compression.decPCK (stream, dim1=None, dim2=None, overflowPix=None, version=None)`

Modified CCP4 pck decompressor used in MAR345 images

Parameters **stream** – string or file

Returns numpy.ndarray (square array)

`fabio.compression.decTY1 (raw_8, raw_16=None, raw_32=None)`

Modified byte offset decompressor used in Oxford Diffraction images

Parameters

- **raw_8** – strings containing raw data with integer 8 bits
- **raw_16** – strings containing raw data with integer 16 bits
- **raw_32** – strings containing raw data with integer 32 bits

Returns numpy.ndarray

`fabio.compression.decZlib (stream)`

Decompress a chunk of data using the zlib algorithm from Python

`fabio.compression.endianness ()`

Return the native endianness of the system

`fabio.compression.md5sum (blob)`

returns the md5sum of an object...

4.30 fabio.converters Module

Converter module. This is for the moment empty (populated only with almost pass through anonymous functions) but aims to be populated with more sophisticated translators ...

`fabio.converters.convert_data (inp, outp, data)`

Return data converted to the output format ... over-simplistic implementation for the moment ... :param inp,outp: input/output format like “cbfimage” :param data(ndarray): the actual dataset to be transformed

`fabio.converters.convert_data_integer (data)`

convert data to integer

`fabio.converters.convert_header (inp, outp, header)`

return header converted to the output format :param inp,outp: input/output format like “cbfimage” :param header(dict):the actual set of headers to be transformed

4.31 fabio.datIO Module

Authors: Henning O. Sorensen & Erik Knudsen Center for Fundamental Research: Metal Structures in Four Dimensions Risoe National Laboratory Frederiksborgvej 399 DK-4000 Roskilde email:erik.knudsen@risoe.dk and Jon Wright, ESRF

class `fabio.datIO.columnfile (data=None, clabels=None, rlabels=None, fname=None)`

Bases: `fabio.datIO.fabiodata`

Concrete fabiodata class

read (*fname, frame=None*)

class `fabio.datIO.fabiodata (data=None, clabels=None, rlabels=None, fname=None)`

Bases: `object`

A common class for dataIO in fable Contains a 2d numpy array for keeping data, and two lists (clabels and rlabels) containing labels for columns and rows respectively

read (*fname=None, frame=None*)

To be overridden by format specific subclasses

4.32 fabio.TiffIO Module

class `fabio.TiffIO.TiffIO (filename, mode=None, cache_length=20, mono_output=False)`

Bases: `object`

getData (*nImage, **kw*)

getImage (*nImage*)

getImageFileDirectories (*fd=None*)

getInfo (*nImage, **kw*)

getNumberOfImages ()

writeImage (*image0, info=None, software=None, date=None*)

4.33 fabio.readbytestream Module

Reads a bytestream

Authors: Jon Wright Henning O. Sorensen & Erik Knudsen ESRF Risoe National Laboratory

`fabio.readbytestream.readbytestream` (*fil*, *offset*, *x*, *y*, *nbytespp*, *datatype*='int', *signed*='n',
swap='n', *typeout*=<type 'numpy.uint16'>)

Reads in a bytestream from a file (which may be a string indicating a filename, or an already opened file (should be "rb")) offset is the position (in bytes) where the pixel data start *nbytespp* = number of bytes per pixel type can be int or float (4 bytes pp) or double (8 bytes pp) signed: normally signed data 'y', but 'n' to try to get back the right numbers when unsigned data are converted to signed (python once had no unsigned numeric types.) swap, normally do not bother, but 'y' to swap bytes typeout is the numpy type to output, normally uint16, but more if overflows occurred *x* and *y* are the pixel dimensions

TODO : Read in regions of interest

PLEASE LEAVE THE STRANGE INTERFACE ALONE - IT IS USEFUL FOR THE BRUKER FORMAT

INDICES AND TABLES

- *genindex*
- *modindex*
- *search*

PYTHON MODULE INDEX

f

- `fabio.__init__`, 15
- `fabio.adscimage`, 22
- `fabio.binaryimage`, 22
- `fabio.bruker100image`, 23
- `fabio.brukerimage`, 23
- `fabio.cbimage`, 24
- `fabio.compression`, 36
- `fabio.converters`, 38
- `fabio.datIO`, 38
- `fabio.dm3image`, 25
- `fabio.edfimage`, 26
- `fabio.fabioimage`, 15
- `fabio.fabioutils`, 17
- `fabio.file_series`, 19
- `fabio.fit2dmaskimage`, 29
- `fabio.fit2dspreadsheetimage`, 29
- `fabio.GEimage`, 29
- `fabio.hdf5image`, 30
- `fabio.HiPiCimage`, 31
- `fabio.kcdimage`, 31
- `fabio.mar345image`, 31
- `fabio.marccdimage`, 32
- `fabio.mrcimage`, 32
- `fabio.openimage`, 22
- `fabio.OXDimage`, 32
- `fabio.pilatusimage`, 33
- `fabio.pixiimage`, 33
- `fabio.pnmimage`, 34
- `fabio.raxisimage`, 34
- `fabio.readbytestream`, 39
- `fabio.TiffIO`, 38
- `fabio.tifimage`, 35
- `fabio.xsdimage`, 35

INDEX

A

acquire() (fabio.fabioutils.DebugSemaphore method), 17
add() (fabio.fabioimage.fabioimage method), 15
adscimage (class in fabio.adscimage), 22
appendFrame() (fabio.edfimage.edfimage method), 27

B

basic_translate() (fabio.brukerimage.brukerimage method), 23
binaryimage (class in fabio.binaryimage), 22
bpp (fabio.edfimage.edfimage attribute), 27
bruker100image (class in fabio.bruker100image), 23
brukerimage (class in fabio.brukerimage), 23
bytecode (fabio.edfimage.edfimage attribute), 27
bytecode (fabio.edfimage.Frame attribute), 26
BZ2File (class in fabio.fabioutils), 17

C

calc_bpp() (fabio.brukerimage.brukerimage method), 23
capsHeader (fabio.edfimage.edfimage attribute), 27
cbfimage (class in fabio.cbfimage), 25
checkData() (fabio.cbfimage.cbfimage static method), 25
checkData() (fabio.fabioimage.fabioimage static method), 15
checkData() (fabio.fit2dmaskimage.fit2dmaskimage static method), 29
checkData() (fabio.kcdimage.kcdimage static method), 31
checkData() (fabio.mar345image.mar345image static method), 31
checkData() (fabio.OXDimage.OXDimage static method), 33
checkData() (fabio.pnmimage.pnmimage static method), 34
checkHeader() (fabio.edfimage.edfimage static method), 27
checkHeader() (fabio.fabioimage.fabioimage static method), 15
CIF (class in fabio.cbfimage), 24
classname (fabio.fabioimage.fabioimage attribute), 15
closed (fabio.fabioutils.GzipFile attribute), 17

columnfile (class in fabio.datIO), 38
compByteOffet_numpy() (in module fabio.compression), 36
compPCK() (in module fabio.compression), 36
compTY1() (in module fabio.compression), 36
construct_filename() (in module fabio.fabioutils), 18
convert() (fabio.fabioimage.fabioimage method), 15
convert_data() (in module fabio.converters), 38
convert_data_integer() (in module fabio.converters), 38
convert_header() (in module fabio.converters), 38
current() (fabio.file_series.file_series method), 19
current() (fabio.file_series.filename_series method), 20
current_image() (fabio.file_series.file_series method), 19
current_image() (fabio.file_series.filename_series method), 20
current_object() (fabio.file_series.file_series method), 19
current_object() (fabio.file_series.filename_series method), 20

D

data (fabio.edfimage.edfimage attribute), 27
data (fabio.edfimage.Frame attribute), 26
DebugSemaphore (class in fabio.fabioutils), 17
decByteOffet_cython() (in module fabio.compression), 36
decByteOffet_numpy() (in module fabio.compression), 36
decByteOffet_python() (in module fabio.compression), 36
decByteOffet_weave() (in module fabio.compression), 37
decBzip2() (in module fabio.compression), 37
decGzip() (in module fabio.compression), 37
decKM4CCD() (in module fabio.compression), 37
deconstruct_filename() (fabio.fabioutils.FilenameObject method), 17
deconstruct_filename() (in module fabio.fabioutils), 18
decPCK() (in module fabio.compression), 37
decTY1() (in module fabio.compression), 37
decZlib() (in module fabio.compression), 37
delCapsHeader() (fabio.edfimage.edfimage method), 27
delData() (fabio.edfimage.edfimage method), 27
deleteFrame() (fabio.edfimage.edfimage method), 27

delHeader() (fabio.edfimage.edfimage method), 27
delHeaderKeys() (fabio.edfimage.edfimage method), 27
demo() (in module fabio.GEimage), 30
demo() (in module fabio.pixiimage), 33
deprecated() (in module fabio.fabioutils), 18
dim1 (fabio.edfimage.edfimage attribute), 27
dim2 (fabio.edfimage.edfimage attribute), 27
dims (fabio.edfimage.edfimage attribute), 27
dm3image (class in fabio.dm3image), 25
do_magic() (in module fabio.openimage), 22

E

edfimage (class in fabio.edfimage), 26
endianness() (in module fabio.compression), 37
estimate_offset_value() (fabio.binaryimage.binaryimage method), 22
exists() (fabio.cbimage.CIF method), 24
existsInLoop() (fabio.cbimage.CIF method), 24
extract_data() (fabio.tifimage.Image_File_Directory_entry method), 35
extract_filename() (in module fabio.fabioutils), 18

F

fabio.__init__ (module), 15
fabio.adscimage (module), 22
fabio.binaryimage (module), 22
fabio.bruker100image (module), 23
fabio.brukerimage (module), 23
fabio.cbimage (module), 24
fabio.compression (module), 36
fabio.converters (module), 38
fabio.datIO (module), 38
fabio.dm3image (module), 25
fabio.edfimage (module), 26
fabio.fabioimage (module), 15
fabio.fabioutils (module), 17
fabio.file_series (module), 19
fabio.fit2dmaskimage (module), 29
fabio.fit2dsheetimage (module), 29
fabio.GEimage (module), 29
fabio.hdf5image (module), 30
fabio.HiPiCimage (module), 31
fabio.kcdimage (module), 31
fabio.mar345image (module), 31
fabio.marccdimage (module), 32
fabio.mrcimage (module), 32
fabio.openimage (module), 22
fabio.OXDimage (module), 32
fabio.pilatusimage (module), 33
fabio.pixiimage (module), 33
fabio.pnmimage (module), 34
fabio.raxisimage (module), 34
fabio.readbytestream (module), 39
fabio.TiffIO (module), 38

fabio.tifimage (module), 35
fabio.xsimage (module), 35
fabiodata (class in fabio.datIO), 38
fabioimage (class in fabio.fabioimage), 15
fastReadData() (fabio.edfimage.edfimage method), 27
fastReadROI() (fabio.edfimage.edfimage method), 27
File (class in fabio.fabioutils), 17
file_series (class in fabio.file_series), 19
filename_series (class in fabio.file_series), 20
FilenameObject (class in fabio.fabioutils), 17
first() (fabio.file_series.file_series method), 19
first_image() (fabio.file_series.file_series method), 19
first_object() (fabio.file_series.file_series method), 19
fit2dmaskimage (class in fabio.fit2dmaskimage), 29
fit2dsheetimage (class in fabio.fit2dsheetimage), 29
Frame (class in fabio.edfimage), 26

G

GEimage (class in fabio.GEimage), 29
gen_header() (fabio.brukerimage.brukerimage method), 23
gen_overflow() (fabio.brukerimage.brukerimage method), 24
getBpp() (fabio.edfimage.edfimage method), 27
getByteCode() (fabio.edfimage.edfimage method), 27
getByteCode() (fabio.edfimage.Frame method), 26
getCapsHeader() (fabio.edfimage.edfimage method), 27
getclassname() (fabio.fabioimage.fabioimage method), 15
getCompressionRatio() (fabio.OXDimage.OXDimage method), 33
getData() (fabio.edfimage.edfimage method), 27
getData() (fabio.edfimage.Frame method), 26
getData() (fabio.TiffIO.TiffIO method), 38
getDim1() (fabio.edfimage.edfimage method), 28
getDim2() (fabio.edfimage.edfimage method), 28
getDims() (fabio.edfimage.edfimage method), 28
getEdfBlock() (fabio.edfimage.Frame method), 26
getframe() (fabio.edfimage.edfimage method), 28
getframe() (fabio.fabioimage.fabioimage method), 15
getframe() (fabio.GEimage.GEimage method), 29
getframe() (fabio.hdf5image.hdf5image method), 30
getframe() (fabio.mrcimage.mrcimage method), 32
getframe() (fabio.pixiimage.pixiimage method), 34
getHeader() (fabio.edfimage.edfimage method), 28
getheader() (fabio.fabioimage.fabioimage method), 15
getHeaderKeys() (fabio.edfimage.edfimage method), 28
getImage() (fabio.TiffIO.TiffIO method), 38
getImageFileDirectories() (fabio.TiffIO.TiffIO method), 38
getInfo() (fabio.TiffIO.TiffIO method), 38
getmax() (fabio.fabioimage.fabioimage method), 15
getmean() (fabio.fabioimage.fabioimage method), 15

getmin() (fabio.fabioimage.fabioimage method), 16
 getNbFrames() (fabio.edfimage.edfimage method), 28
 getnum() (in module fabio.fabioutils), 18
 getNumberOfImages() (fabio.TiffIO.TiffIO method), 38
 getSize() (fabio.fabioutils.BZ2File method), 17
 getSize() (fabio.fabioutils.File method), 17
 getSize() (fabio.fabioutils.GzipFile method), 17
 getSize() (fabio.fabioutils.StringIO method), 18
 getSize() (fabio.OXDImage.Section method), 33
 getstddev() (fabio.fabioimage.fabioimage method), 16
 GzipFile (class in fabio.fabioutils), 17

H

hdf5image (class in fabio.hdf5image), 30
 HDF5location (class in fabio.hdf5image), 30
 header (fabio.edfimage.edfimage attribute), 28
 header_keys (fabio.edfimage.edfimage attribute), 28
 HiPiCimage (class in fabio.HiPiCimage), 31

I

Image_File_Directory (class in fabio.tifimage), 35
 Image_File_Directory_entry (class in fabio.tifimage), 35
 integrate_area() (fabio.fabioimage.fabioimage method), 16
 interpret_header() (in module fabio.marccdimage), 32
 isAscii() (fabio.cbimage.CIF static method), 24
 isAscii() (in module fabio.fabioutils), 18

J

jump() (fabio.file_series.file_series method), 19
 jump() (fabio.file_series.filename_series method), 20
 jump_filename() (in module fabio.fabioutils), 18
 jump_image() (fabio.file_series.file_series method), 19
 jump_image() (fabio.file_series.filename_series method), 20
 jump_object() (fabio.file_series.file_series method), 20
 jump_object() (fabio.file_series.filename_series method), 20

K

kcdimage (class in fabio.kcdimage), 31

L

last() (fabio.file_series.file_series method), 20
 last_image() (fabio.file_series.file_series method), 20
 last_object() (fabio.file_series.file_series method), 20
 len() (fabio.file_series.file_series method), 20
 load() (fabio.fabioimage.fabioimage method), 16
 loadCHIPLOT() (fabio.cbimage.CIF method), 24
 loadCIF() (fabio.cbimage.CIF method), 24
 LoopHasKey() (fabio.cbimage.CIF static method), 24

M

make_format() (in module fabio.marccdimage), 32

make_slice() (fabio.fabioimage.fabioimage method), 16
 mar345image (class in fabio.mar345image), 31
 marccdimage (class in fabio.marccdimage), 32
 md5sum() (in module fabio.compression), 37
 mrcimage (class in fabio.mrcimage), 32

N

nb_overflow_pixels() (fabio.mar345image.mar345image method), 31
 new_file_series() (in module fabio.file_series), 21
 new_file_series0() (in module fabio.file_series), 21
 next() (fabio.edfimage.edfimage method), 28
 next() (fabio.fabioimage.fabioimage method), 16
 next() (fabio.file_series.file_series method), 20
 next() (fabio.file_series.filename_series method), 21
 next() (fabio.GEimage.GEimage method), 29
 next() (fabio.hdf5image.hdf5image method), 30
 next() (fabio.mrcimage.mrcimage method), 32
 next() (fabio.pxiimage.pxiimage method), 34
 next_filename() (in module fabio.fabioutils), 18
 next_image() (fabio.file_series.file_series method), 20
 next_image() (fabio.file_series.filename_series method), 21
 next_object() (fabio.file_series.file_series method), 20
 next_object() (fabio.file_series.filename_series method), 21
 nframes (fabio.edfimage.edfimage attribute), 28
 nice_int() (in module fabio.fabioutils), 18
 numbered_file_series (class in fabio.file_series), 21
 numstem() (in module fabio.fabioutils), 18

O

openheader() (in module fabio.openimage), 22
 openimage() (in module fabio.openimage), 22
 OXDImage (class in fabio.OXDImage), 33

P

P1dec() (fabio.pnmimage.pnmimage method), 34
 P2dec() (fabio.pnmimage.pnmimage method), 34
 P3dec() (fabio.pnmimage.pnmimage static method), 34
 P4dec() (fabio.pnmimage.pnmimage static method), 34
 P5dec() (fabio.pnmimage.pnmimage method), 34
 P6dec() (fabio.pnmimage.pnmimage static method), 34
 P7dec() (fabio.pnmimage.pnmimage static method), 34
 pad() (in module fabio.fabioutils), 18
 parse() (fabio.hdf5image.HDF5location method), 30
 parseheader() (fabio.edfimage.Frame method), 26
 pilatusimage (class in fabio.pilatusimage), 33
 pxiimage (class in fabio.pxiimage), 33
 pnmimage (class in fabio.pnmimage), 34
 pop() (fabio.cbimage.CIF method), 25
 popitem() (fabio.cbimage.CIF method), 25
 prev_image() (fabio.file_series.filename_series method), 21

previous() (fabio.edfimage.edfimage method), 28
 previous() (fabio.fabioimage.fabioimage method), 16
 previous() (fabio.file_series.file_series method), 20
 previous() (fabio.file_series.filename_series method), 21
 previous() (fabio.GEimage.GEimage method), 30
 previous() (fabio.hdf5image.hdf5image method), 30
 previous() (fabio.mrcimage.mrcimage method), 32
 previous() (fabio.pixiimage.pixiimage method), 34
 previous_filename() (in module fabio.fabioutils), 19
 previous_image() (fabio.file_series.file_series method), 20
 previous_object() (fabio.file_series.file_series method), 20
 previous_object() (fabio.file_series.filename_series method), 21

R

raxisimage (class in fabio.raxisimage), 34
 read() (fabio.adscimage.adscimage method), 22
 read() (fabio.binaryimage.binaryimage method), 23
 read() (fabio.bruker100image.bruker100image method), 23
 read() (fabio.brukerimage.brukerimage method), 24
 read() (fabio.cbimage.cbimage method), 25
 read() (fabio.datIO.columnfile method), 38
 read() (fabio.datIO.fabiodata method), 38
 read() (fabio.dm3image.dm3image method), 25
 read() (fabio.edfimage.edfimage method), 28
 read() (fabio.fabioimage.fabioimage method), 16
 read() (fabio.fit2dmaskimage.fit2dmaskimage method), 29
 read() (fabio.fit2dspreadsheetimage.fit2dspreadsheetimage method), 29
 read() (fabio.GEimage.GEimage method), 30
 read() (fabio.hdf5image.hdf5image method), 30
 read() (fabio.HiPiCimage.HiPiCimage method), 31
 read() (fabio.kcdimage.kcdimage method), 31
 read() (fabio.mar345image.mar345image method), 31
 read() (fabio.mrcimage.mrcimage method), 32
 read() (fabio.OXDimage.OXDimage method), 33
 read() (fabio.pixiimage.pixiimage method), 34
 read() (fabio.pnmimage.pnmimage method), 34
 read() (fabio.raxisimage.raxisimage method), 35
 read() (fabio.tifimage.tifimage method), 35
 read() (fabio.xsdimage.xsdimage method), 36
 read_data() (fabio.dm3image.dm3image method), 25
 read_tag_entry() (fabio.dm3image.dm3image method), 25
 read_tag_group() (fabio.dm3image.dm3image method), 25
 read_tag_type() (fabio.dm3image.dm3image method), 25
 readbytes() (fabio.dm3image.dm3image method), 26
 readbytestream() (in module fabio.readbytestream), 39
 readCIF() (fabio.cbimage.CIF method), 25

readheader() (fabio.fabioimage.fabioimage method), 16
 readROI() (fabio.fabioimage.fabioimage method), 16
 rebin() (fabio.fabioimage.fabioimage method), 16
 release() (fabio.fabioutils.DebugSemaphore method), 17
 resetvals() (fabio.fabioimage.fabioimage method), 16
 rigakuKeys() (fabio.raxisimage.raxisimage method), 35

S

save() (fabio.fabioimage.fabioimage method), 16
 saveCIF() (fabio.cbimage.CIF method), 25
 Section (class in fabio.OXDimage), 33
 seek() (fabio.fabioutils.GzipFile method), 17
 set_index() (fabio.hdf5image.HDF5location method), 30
 setBpp() (fabio.edfimage.edfimage method), 28
 setByteCode() (fabio.edfimage.edfimage method), 28
 setByteCode() (fabio.edfimage.Frame method), 26
 setCapsHeader() (fabio.edfimage.edfimage method), 28
 setData() (fabio.edfimage.edfimage method), 28
 setData() (fabio.edfimage.Frame method), 26
 setData() (fabio.OXDimage.Section method), 33
 setDim1() (fabio.edfimage.edfimage method), 28
 setDim2() (fabio.edfimage.edfimage method), 28
 setHeader() (fabio.edfimage.edfimage method), 28
 setHeaderKeys() (fabio.edfimage.edfimage method), 28
 setNbFrames() (fabio.edfimage.edfimage method), 28
 setSize() (fabio.fabioutils.BZ2File method), 17
 setSize() (fabio.fabioutils.File method), 17
 setSize() (fabio.fabioutils.GzipFile method), 18
 setSize() (fabio.fabioutils.StringIO method), 18
 size (fabio.fabioutils.BZ2File attribute), 17
 size (fabio.fabioutils.File attribute), 17
 size (fabio.fabioutils.GzipFile attribute), 18
 size (fabio.fabioutils.StringIO attribute), 18
 str() (fabio.fabioutils.FilenameObject method), 17
 StringIO (class in fabio.fabioutils), 18
 swap_needed() (fabio.binaryimage.binaryimage static method), 23
 swap_needed() (fabio.edfimage.edfimage method), 28
 swap_needed() (fabio.edfimage.Frame method), 26
 swap_needed() (fabio.raxisimage.raxisimage method), 35

T

test() (in module fabio.adscimage), 22
 test() (in module fabio.brukerimage), 24
 test() (in module fabio.fabioimage), 16
 Tiff_header (class in fabio.tifimage), 35
 TiffIO (class in fabio.TiffIO), 38
 tifimage (class in fabio.tifimage), 35
 to_url() (fabio.hdf5image.HDF5location method), 30
 toAscii() (in module fabio.fabioutils), 19
 toPIL16() (fabio.bruker100image.bruker100image method), 23
 toPIL16() (fabio.fabioimage.fabioimage method), 16
 toString() (fabio.cbimage.CIF method), 25

tostring() (fabio.fabioutils.FilenameObject method), 17

U

UnknownCompressedFile (class in fabio.fabioutils), 18

unpack() (fabio.edfimage.edfimage method), 29

unpack() (fabio.tifimage.Image_File_Directory method),
35

unpack() (fabio.tifimage.Image_File_Directory_entry
method), 35

update_header() (fabio.fabioimage.fabioimage method),
16

W

write() (fabio.adscimage.adscimage method), 22

write() (fabio.binaryimage.binaryimage method), 23

write() (fabio.bruckerimage.bruckerimage method), 24

write() (fabio.cbimage.cbimage method), 25

write() (fabio.edfimage.edfimage method), 29

write() (fabio.fabioimage.fabioimage method), 16

write() (fabio.fit2dmaskimage.fit2dmaskimage method),
29

write() (fabio.GEimage.GEimage method), 30

write() (fabio.hdf5image.hdf5image method), 30

write() (fabio.mar345image.mar345image method), 31

write() (fabio.mrcimage.mrcimage method), 32

write() (fabio.OXDimage.OXDimage method), 33

write() (fabio.pixiimage.pixiimage method), 34

write() (fabio.pnmimage.pnmimage method), 34

write() (fabio.tifimage.tifimage method), 35

writeImage() (fabio.TiffIO.TiffIO method), 38

X

xsdimage (class in fabio.xsdimage), 35